

Technical Bulletin

Chemical Compatibility

| Chemical | Container Materials | | | | | Closure Liner Materials | | | | | Closure Materials | | | | Septa and Stopper Materials | | | | | | |
|--------------------------------|---------------------|------|------|-----|----|-------------------------|------|------------|-----|----------|-------------------|-----|----------|----|-----------------------------|--------------|------|----------------|----------|------|-----|
| | Glass | HDPE | LDPE | PET | PP | Al Foil | LDPE | Poly-Vinyl | SBR | Silicone | PTFE | PBT | Phenolic | PP | Urea | Butyl Rubber | EPDM | Natural Rubber | Silicone | PTFE | FKM |
| Acetic Acid, Gglacial | A | A | B | A | A | A | B | B | C | B | A | C | A | A | D | B | B | B | B | A | D |
| Acetone | A | D | D | C | B | A | D | D | D | B | A | D | A | B | A | B | A | B | B | A | D |
| Acetonitrile | A | A | A | B | A | A | A | D | B | D | A | — | A | A | — | D | B | D | D | A | D |
| Acrylonitrile | A | A | A | B | B | B | A | D | C | D | A | — | D | B | — | D | A | D | D | A | D |
| Ammonium Sulfide | A | A | A | — | A | D | A | A | B | A | A | — | A | A | C | A | A | A | A | A | C |
| Benzene | A | D | D | C | D | B | D | D | D | D | A | A | A | D | A | D | D | D | D | A | A |
| Bleach | A | A | B | C | B | D | B | A | D | B | A | C | D | B | — | A | A | D | B | A | A |
| Boric Acid | A | A | A | A | A | D | A | A | A | A | A | A | B | A | — | A | A | A | A | A | A |
| Carbonic Acid | A | A | A | — | A | B | A | A | B | A | A | — | — | A | — | A | A | A | A | A | A |
| Chlorobenzene | A | C | D | B | C | A | D | D | D | D | A | B | A | C | B | D | D | D | D | A | A |
| Chloroform | A | C | C | D | D | A | C | D | D | D | A | D | A | D | A | D | D | D | D | A | A |
| Dichloromethane (DCM) | A | C | D | D | C | D | D | D | D | D | A | D | C | C | B | D | C | D | D | A | B |
| Diethylamine | A | C | D | — | B | A | D | D | B | B | A | — | — | B | — | B | B | B | B | A | C |
| Dimethyl Formamide (DMF) | A | A | A | B | A | A | A | D | D | B | A | — | A | A | — | D | B | D | B | A | D |
| Dimethyl Sulfoxide (DMSO) | A | A | A | B | A | A | A | D | D | D | A | — | — | A | — | D | A | D | D | A | D |
| Dioxane | A | B | B | A | D | D | B | D | D | D | A | B | A | D | — | B | B | D | D | A | D |
| Ether | A | C | D | A | D | B | D | D | D | D | A | — | B | D | B | D | C | D | D | A | C |
| Ethyl Acetate | A | B | B | B | C | B | B | D | D | C | A | C | A | C | B | C | C | D | C | A | D |
| Ethyl Alcohol | A | A | A | A | A | B | A | B | A | B | A | A | B | A | A | A | A | A | B | A | A |
| Ethylene Glycol | A | A | A | A | A | B | A | A | A | A | A | A | B | A | B | A | A | A | A | A | A |
| Formaldehyde | A | A | A | B | A | A | A | C | B | B | A | — | B | A | A | A | A | C | B | A | C |
| Formic Acid 50% | A | A | B | — | A | C | B | B | B | C | A | A | C | A | D | A | A | B | C | A | C |
| Gasoline | A | C | D | B | C | A | D | D | D | D | A | A | B | C | A | D | D | D | D | A | A |
| Glycerine | A | A | A | — | A | A | A | C | A | B | A | — | A | A | — | A | A | A | B | A | A |
| Heptane | A | C | D | B | C | A | D | C | D | D | A | A | A | C | A | D | D | D | D | A | A |
| Hexane | A | B | D | C | B | A | D | D | D | D | A | A | B | B | — | D | D | D | D | A | A |
| Hydrochloric Acid (HCL) 50% | A | A | A | B | A | D | A | B | D | D | A | C | A | A | D | A | D | B | D | A | A |
| Hydrofluoric Acid (HF) 50% | D | A | A | C | A | D | A | C | D | D | A | C | D | A | D | C | D | C | D | A | A |
| Hydrogen Peroxide 50% | B | A | A | B | A | A | A | C | C | B | A | B | D | A | D | B | B | B | B | A | A |
| Iodine | A | C | D | A | C | A | D | C | B | A | A | — | — | C | — | B | B | D | A | A | A |
| Isopropyl Alcohol | A | A | A | A | A | A | A | B | B | A | A | A | A | A | — | A | A | A | A | A | A |
| Methyl Alcohol | A | A | A | B | A | A | A | C | A | A | A | B | B | A | A | A | A | A | A | A | D |
| Methyl Ethyl Ketone (MEK) | A | D | D | B | B | A | D | D | D | D | A | C | A | B | — | A | B | D | D | A | D |
| Methylene Chloride | A | C | D | D | C | D | D | D | D | D | A | D | C | C | B | D | D | D | D | A | B |
| Nitric Acid 50% | A | C | B | C | C | D | B | B | D | D | A | C | B | C | D | C | D | C | D | A | B |
| Pentane | A | C | C | — | D | A | C | D | D | D | A | — | — | D | — | D | D | D | D | A | A |
| Perchloric Acid 50% | B | B | B | B | B | D | B | D | D | D | B | — | — | B | — | B | B | D | D | B | A |
| Phenol 50% | A | D | D | D | D | A | D | C | D | D | A | D | A | D | — | D | C | D | D | A | A |
| Phosphoric Acid 50% | A | A | A | B | A | B | A | B | D | D | A | B | B | A | D | B | B | D | D | A | A |
| Picric Acid | A | D | D | B | D | A | D | D | B | D | A | — | A | D | D | B | B | B | D | A | A |
| Potassium Hydroxide | D | A | A | D | A | D | A | A | B | C | A | C | D | A | — | A | A | B | C | A | B |
| Sodium Hydroxide 50% | D | A | B | D | A | D | B | C | A | B | A | C | D | A | C | A | A | A | B | A | B |
| Sodium Peroxide | A | B | B | — | B | C | B | A | B | D | A | — | B | B | D | A | A | B | D | A | A |
| Sodium Thiosulfate | A | A | A | B | A | A | A | A | B | A | A | — | A | A | B | A | A | B | A | A | A |
| Sulfuric Acid 50% | A | A | A | B | B | C | A | C | D | D | A | B | C | B | D | D | C | D | D | A | A |
| Tetrahydrofuran (THF) | A | C | C | A | B | A | C | D | D | D | A | D | A | B | — | C | C | D | D | A | D |
| Toluene | A | C | C | C | C | A | C | C | D | D | A | D | A | C | — | D | D | D | D | A | B |
| Trifluoroacetic acid (TFA) 50% | A | A | A | B | A | — | A | A | — | D | — | — | — | A | — | — | A | — | D | — | — |
| Vegetable Oil | A | B | B | A | C | A | B | A | D | A | A | A | A | A | A | C | C | D | A | A | A |
| Xylene | A | C | D | C | D | A | D | D | D | D | A | C | A | D | B | D | D | D | D | A | A |

(Tests conducted at room temp) **A** - Resistant **B** - Limited Resistance **C** - Poor Resistance **D** - Not Resistant **—** - Unknown
 Al Foil ... aluminum foil EPDM ... ethylene propylene diene monomer FMK ... fluorocastomer HDPE ... high density polyethylene LDPE ... low density polyethylene PBT ... polybutylene terephthalate
 PET ... polyethylene terephthalate PP ... polypropylene PTFE ... polytetrafluoroethylene SBR ... styrene butadiene rubber Urea ... urea

Although the information in this chart was acquired from reputable sources, it should only be used as a guide in selecting a container and closure system. Because so many factors can affect the chemical resistance of a material, in-house testing under actual conditions should be performed. WHEATON accepts no responsibility for the accuracy of this data or for any consequences resulting from its use.