

C&G Containers Scientific

"The Total Package!"



EPA UCMR5 Compliant Certified PFAS Containers

Minimize contamination risks & maximize accuracy in your test results.

METHOD 1633

250cc Natural HDPE Oblong LEAKPROOF 45PP cap
Part #LPV008200650

METHOD 537.1

250cc Natural HDPE Oblong LEAKPROOF 45PP cap
TRIS/TRIS HCL (1.25gms)
Part #LPV008249600

METHOD 533

250cc Natural HDPE Oblong LEAKPROOF 45PP cap
Ammonia Acetate (250mg)
Part #LPV008294600

(800) 369-5524
sales@foxscientific.com
www.foxscientific.com


Fox
Scientific, Inc.
www.foxscientific.com • sales@foxscientific.com
(800) 369-5524