



**WHITE DIAMOND SERIES
SMALL VOLUME STORAGE
REFRIGERATORS &
FREEZERS**



8221 East FM 917 • Alvarado, TX 76009 • Phone:(800)369-5524 • Fax:(817)783-3571
www.foxscientific.com • sales@foxscientific.com



Product Portfolio

White Diamond Series small volume storage product offering includes 20+ refrigerators and 15+ freezers between 1 and 10.5 ft³. Freezers are available with either manual, auto, or controlled auto defrost options.

Environmentally Friendly

White Diamond Series small volume storage product offering utilizes natural hydrocarbon refrigerants that vastly reduce global-warming potential while meeting EPA/SNAP mandates.

Energy Star Certified

Corepoint™ Scientific offers 20+ undercounter and compact refrigerators and freezers that meet and/or exceed Energy Star Certification standards.

UL/UC-L Listed

White Diamond Series small volume storage product offering has been tested by the Underwriters Laboratories and have met and exceeded all UL/UC-L listing requirements.

WHY CHOOSE US ?

Corepoint™ Scientific's White Diamond Series undercounter, compact, and dual temperature combination refrigerators and freezers are designed to meet the demanding requirements for scientific and laboratory research. Our advanced engineered design incorporates the latest in cabinet refrigeration, temperature control and monitoring features.

The White Diamond Series provides energy efficient, convenient, safe and reliable performance for optimal storage temperature environments necessary for a wide range of life science, pharmacy, biological, medical, clinical, vaccine and industrial applications.

- *Digital temperature displays*
- *Integrated alarm systems*
- *Interior lighting (select models only)*
- *Probe access*
- *Freestanding or Built-in Configurations Available*



Freestanding Refrigerator	Ft'	Set Point	Defrost	Door	# of Shelves	Dimensions (w" x d" x h")
*PR011WWW/0	1	4°C	Cycle	1 Solid	2 Adj.	17 ¼ x 19 ¼ x 21 ¼
*PR011WWG/0	1	4°C	Cycle	1 Glass	2 Adj.	17 ¼ x 19 ¼ x 21 ¼
PR031WWW/0	2.5	4°C	Cycle	1 Solid	4 Adj.	17 ¼ x 19 ½ x 28 ¾
PR031WWG/0	2.5	4°C	Cycle	1 Glass	4 Adj.	17 ¼ x 19 ½ x 28 ¾
PR061WWW/0	5.2	4°C	Cycle	1 Solid	3 Adj.	23 ¾ x 24 x 32 ½
PR061WWG/0	5.2	4°C	Cycle	1 Glass	3 Adj.	23 ¾ x 24 x 32 ½
PR101WWW/0	10.5	4°C	Cycle	1 Solid	6 Adj.	20 ¾ x 26 ½ x 59
PR101WWG/0	10.5	4°C	Cycle	1 Glass	6 Adj.	20 ¾ x 26 ½ x 59
Freestanding Freezer	Ft'	Set Point	Defrost	Door	# of Shelves	Dimensions (w" x d" x h")
PF021WWW/0CAD	1	-20°C	Controlled Auto	1 Solid	2 Adj.	23 ¾ x 24 x 20 ¾
*PF021WWW/0A	1.3	-20°C	Auto	1 Solid	2 Adj.	23 ¾ x 24 x 20 ¾
*PF021WWW/0M	1.7	-20°C	Manual	1 Solid	2 Adj.	23 ¾ x 24 x 20 ¾
PF051WWW/0M20	4	-20°C	Manual	1 Solid	2 Fixed	23 ¾ x 24 ¼ x 32 ¾
PF051WWW/0M30	4	-30°C	Manual	1 Solid	3 Adj.	23 ¾ x 24 x 32 ¾
PF051WWW/0M40	4	-40°C	Manual	1 Solid	3 Adj.	23 ¾ x 24 x 32 ¾
Freestanding Combo Model	Ft' (Ref/Frz)	Set Point	Defrost (Ref/Frz)	Door (Ref/Frz)	# of Shelves (Ref/Frz)	Dimensions (w" x d" x h")
PRF072WWG/0CAD	5.2/1	4°C/-20°C	Cycle/Controlled Auto	1 Glass/1 Solid	3 Adj./2 Adj.	23 ¾ x 24 ¼ x 53 ¾
PRF072WWW/0CAD	5.2/1	4°C/-20°C	Cycle/Controlled Auto	2 Solid	3 Adj./2 Adj.	23 ¾ x 24 x 53 ¾
PRF072WWG/0A	5.2/1.3	4°C/-20°C	Cycle/Auto	1 Glass/1 Solid	3 Adj./2 Adj.	23 ¾ x 24 x 53 ¾
PRF072WWW/0A	5.2/1.3	4°C/-20°C	Cycle/Auto	2 Solid	3 Adj./2 Adj.	23 ¾ x 24 x 53 ¾
PRF072WWG/0M	5.2/1.7	4°C/-20°C	Cycle/Manual	1 Glass/1 Solid	3 Adj./2 Adj.	23 ¾ x 24 x 53 ¾
PRF072WWW/0M	5.2/1.7	4°C/-20°C	Cycle/Manual	2 Solid	3 Adj./2 Adj.	23 ¾ x 24 x 53 ¾
PRF102WWG/0	5.2/4	4°C/-20°C	Cycle/Manual	1 Glass/1 Solid	3 Adj./2 Fixed	23 ¾ x 24 x 65
PRF102WWW/0	5.2/4	4°C/-20°C	Cycle/Manual	2 Solid	3 Adj./2 Fixed	23 ¾ x 24 x 65
PRF102WWG/030	5.2/4	4°C/-30°C	Cycle/Manual	1 Glass/1 Solid	3 Adj./3 Adj.	23 ¾ x 24 x 65
PRF102WWW/030	5.2/4	4°C/-30°C	Cycle/Manual	2 Solid	3 Adj./3 Adj.	23 ¾ x 24 x 65
PRF102WWG/040	5.2/4	4°C/-40°C	Cycle/Manual	1 Glass/1 Solid	3 Adj./3 Adj.	23 ¾ x 24 x 65
PRF102WWW/040	5.2/4	4°C/-40°C	Cycle/Manual	2 Solid	3 Adj./3 Adj.	23 ¾ x 24 x 65
PRF122WWG/0CAD	10.5/1	4°C/-20°C	Cycle/Controlled Auto	1 Glass/1 Solid	6 Adj./2 Adj.	23 ¾ x 26 ½ x 80 ¾
PRF122WWW/0CAD	10.5/1	4°C/-20°C	Cycle/Controlled Auto	2 Solid	6 Adj./2 Adj.	23 ¾ x 26 ½ x 80 ¾
PRF122WWG/0A	10.5/1.3	4°C/-20°C	Cycle/Auto	1 Glass/1 Solid	6 Adj./2 Adj.	23 ¾ x 26 ½ x 80 ¾
PRF122WWW/0A	10.5/1.3	4°C/-20°C	Cycle/Auto	2 Solid	6 Adj./2 Adj.	23 ¾ x 26 ½ x 80 ¾
PRF122WWG/0M	10.5/1.7	4°C/-20°C	Cycle/Manual	1 Glass/1 Solid	6 Adj./2 Adj.	23 ¾ x 26 ½ x 80 ¾
PRF122WWW/0M	10.5/1.7	4°C/-20°C	Cycle/Manual	2 Solid	6 Adj./2 Adj.	23 ¾ x 26 ½ x 80 ¾
Built-in Refrigerator	Ft'	Set Point	Defrost	Door	# of Shelves	Dimensions (w" x d" x h")
*PR031WWW/0FB	2.5	4°C	Cycle	1 Solid	4 Adj.	17 ¾ x 22 ¾ x 30 ¾
*PR051SSS/0	4.6	4°C	Cycle	1 Solid	2 Adj./1 Fixed	23 ¾ x 24 ½ x 33 ¾
*PR051WWW/0	4.6	4°C	Cycle	1 Solid	2 Adj./1 Fixed	23 ¾ x 24 ½ x 33 ¾
*PR051WWG/0	4.6	4°C	Cycle	1 Glass	2 Adj./1 Fixed	23 ¾ x 25 ½ x 33 ¾
*PR051WWWADA/0	4.6	4°C	Cycle	1 Solid	2 Adj./1 Fixed	23 ¾ x 24 ½ x 31 ½ ¹⁶
*PR051WWGADA/0	4.6	4°C	Cycle	1 Glass	2 Adj./1 Fixed	23 ¾ x 25 ½ x 31 ½ ¹⁶
Built-in Freezer	Ft'	Set Point	Defrost	Door	# of Shelves	Dimensions (w" x d" x h")
*PF051WWW/0M	4.2	-20°C	Manual	1 Solid	2 Fixed	23 ¾ x 24 ½ x 33 ¾
*PF051WWWADA/0M	4.2	-20°C	Manual	1 Solid	2 Fixed	23 ¾ x 24 ½ x 31 ½ ¹⁶
*PF051SSS/0M	4.2	-20°C	Manual	1 Solid	2 Fixed	23 ¾ x 24 ½ x 33 ¾
*PF051WWW/0A	4.2	-20°C	Auto	1 Solid	3 Adj.	23 ¾ x 24 ½ x 33 ¾
Built-In Combo Model	Ft' (Ref/Frz)	Set Point (Ref/Frz)	Defrost (Ref/Frz)	Door (Ref/Frz)	# of Shelves (Ref/Frz)	Dimensions (w" x d" x h")
*PRF092WWG/0	4.6/4.2	4°C/-20°C	Cycle/Manual	1 Glass/1 Solid	2 Adj., 1 Fixed/2 Fixed	23 ¾ x 25 ½ x 66 ¾
*PRF092WWW/0	4.6/4.2	4°C/-20°C	Cycle/Manual	2 Solid	2 Adj., 1 Fixed/2 Fixed	23 ¾ x 24 ¾ x 66 ¾
*PRF092SSS/0	4.5/4.2	4°C/-20°C	Cycle/Manual	2 Solid	2 Adj., 1 Fixed/2 Fixed	23 ¾ x 24 ½ x 66 ¾

*Left hinge model available