

Key Products for Environmental Testing

HPLC Columns



79181 - PRP-X100 Anion Exchange HPLC Column, 4.6 x 250 mm, 5 μ m, PEEK



79443 - PRP-1 Reversed-Phase HPLC Column, 4.1 x 50 mm, 5 μ m, Stainless Steel



79676 - PRP-C18 Reversed-Phase HPLC Column, 4.6 x 150 mm, 5 μ m, Stainless Steel

HPLC Autosampler Syringes

For: Agilent 1290 Infinity LC Injector HTC HTS, Waters 2777 Sample Manager, UltiMate 3000, PAL LC HTC and HTS:

203073: 10 μ L Model 701 FN CTC SYR (6.6 mm) C-Line, Cemented NDL, 22s gauge, 2 in, point style 3



For: PAL LC HTC and HTS Syringes, Ulti-Mate 3000 Syringes, Agilent 1290 Infinity LC Injector HTC HTS Syringes, Waters 2777 Sample Manager

203235 - 100 μ L, Model 1710 FN CTC SYR (6.6 mm), C-Line, Cemented NDL, 22 gauge, point style 3.

Microlab 600 Diluters and Dispensers

ML615-CNT

Dual Syringe Continuous Dispenser with Basic Controller and Concorde Hand Probe



ML625-DTHP

Dual Syringe Diluter with Advanced Controller and Disposable Tip Hand Probe



Manual Syringes



86336 - 50 mL Gastight Model 1050 Sample Lock (SL), Large Removable Needle, 22 gauge, 2 in, point style 2



86020 - 100 mL Gastight Model 1100, PTFE Luer Lock (TLL), Needle Sold Separately

Pipettes



55019-30 - 0.2 to 2 μ L adjustable volume pipette



55019-40 - 100 to 1000 μ L adjustable volume pipette

GC Autosampler Syringes



For: Agilent GC Autosamplers 7693A, 7673, 7683, and 6850 ALS
80390 - 10 μ L Microliter Agilent Autosampler Syringe Model 701 N, Cemented Needle, 23s gauge, 1.71 in, point style AS (6 pack).



For: Syringe for Shimadzu AOC-14 AOC-17 and AOC-20
93898-01 - 10 μ L Microliter Shimadzu Syringe Model 701 N, Cemented Needle, 23s gauge, 1.69 in, point style AS.

