


Granite Platforms with Vibration Isolators 2 sizes now available for use with Sartorius Balances.

Sartorius now offers granite platforms with specially designed silicone vibration isolators. This new accessory is ideal for use with analytical, semi-micro, micro and ultra-micro balances. These platforms will help eliminate vibration and enhance balance response time.

Standard Model for use with most Sartorius Balances

The standard model consists of a 13" x 15" granite platform. The isolators are a high-grade silicone connected to the platform and the base. These specially designed silicone isolators filter out vibrations to improve balance performance.

13" x 15" Anti-Vibration Platform Sartorius Catalog # 21201315



NEW Extended Model for use with Sartorius ME Genius Series Balances

The extended model consists of a 16" x 21" granite platform. The isolators are a high-grade silicone connected to the platform and the base. Designed for use with larger balances, such as the Sartorius GENIUS Series balances.

16" x 21" Anti-Vibration Platform Sartorius Catalog # 24201621

Fox Scientific, Inc

8221 East FM 917
Alvarado, TX 76009
Toll Free: (800) 369-5524
Phone (817)783-7500 Fax: (817) 783-3571
sales@foxscientific.com
www.foxscientific.com