

Tough-Tags® Tube Labels

Pre-cut labels for identification

- Adheres to most plastics, glass, and metals without cracking, peeling or degrading
- Flexible polyolefin material and acrylic adhesive



Tough-Tags® on a Roll




Ideal Use	0.2 mL Tubes	0.5 mL Tubes	1.5/2.0 mL Tubes	General Purpose	Microplates
Label Diameter	0.81 x 0.28 in/21 x 7 mm	0.94 x 0.50 in/24 x 13 mm	1.28 x 0.50 in/33 x 13 mm	1.50 x 0.75 in/38 x 19 mm	1.50 x 0.28 in/38 x 6 mm
UOM	1,500/roll	1,000/roll	1,000/roll	1,000/roll	1,000/roll
Item No.					
White	HEATT-TNY	HEATTSW-1000	HEATTLW-1000	HEATTGP-1000	HEATT-SWALL
Blue	-	HEATTSB-1000	HEATTLB-1000	-	-
Green	-	HEATTSG-1000	HEATTLG-1000	-	-
Orange	-	HEATTSO-1000	HEATTL0-1000	-	-
Red	-	HEATTSR-1000	HEATTLR-1000	-	-
Yellow	-	HEATTSY-1000	HEATTLY-1000	-	-
Clear	-	HEATTSC-1000	HEATTLc-1000	-	-

- Label sheets are top coated and heat-resistant to prevent jams in laser printers
- Printer formatting instruction included



Tough-Tags®
is a registered trademark
of Diversified Biotech, Inc

Laser Tough-Tags® on Sheets formatted for laser printers

Ideal Use	0.5 mL Tubes	1.5/2.0 mL Tubes	General Purpose	Microscope Slides
				
Label Dimensions	0.94 x 0.50 in/24 x 13 mm	1.28 x 0.50 in/33 x 13 mm	1.50 x 0.75 in/38 x 19 mm	0.94 x 0.94 in/22 x 22 mm
Labels/Sheet	119	85	60	96
Labels/pk	2,975	2,125	1,500	2,400
UOM	25/pk	25/pk	25/pk	25/pk
Item No.				
White	HEATTSW-2240	HEATTLW-2016	HEATTLW-1050	HEAMISL-1000
Blue	HEATTSW-2240-B	HEATTLW-2016-B	HEATTLW-1050-B	HEAMISL-1000-B
Green	HEATTSW-2240-G	HEATTLW-2016-G	HEATTLW-1050-G	HEAMISL-1000-G
Orange	HEATTSW-2240-O	HEATTLW-2016-O	HEATTLW-1050-O	HEAMISL-1000-O
Red	HEATTSW-2240-R	HEATTLW-2016-R	HEATTLW-1050-R	HEAMISL-1000-R
Yellow	HEATTSW-2240-Y	HEATTLW-2016-Y	HEATTLW-1050-Y	HEAMISL-1000-Y
Rainbow*	HEARNBW-3100	HEARNBW-3000	HEARNBW-3200	HEARNBW-3300

*Rainbow 4 sheets of each Blue, Green, Orange, Red, Yellow

Label Station

Organize your labels for easy accessibility


Portable label station holds up to six standard size boxes of Tough-Tags®, Cryo-Babies® and Tough-Spots®

Just drop the box in the label station and feed the paper through the slot

- Acrylic station includes a writing surface, a cutting edge for easy dispensing and non-skid rubber feet
- Counter-sunk holes provided for wall mounting

Item No.	Color	L x W x H in	L x W x H mm
HS23460	Blue	11 x 4.7 x 3.5	28 x 12 x 8.9

Tough-Tags®, Cryo-Babies® and Tough-Spots®



Labels shown not included,
order separately