



Fox Chemicals • Laboratory Supplies • Equipment
Scientific, Inc. est. 1988

OAKTON[®]
I N S T R U M E N T S

**Buy 1 WD-35634-35
at 15% off and
get the WD-35653-04
FREE!**



The Oakton pH pocket tester stands up to everyday use, harsh environments and the test of time. The tester can even stand upright when the cap I used as the stabilizing base.

Double-junction pH electrode sensor is dependable and long lasting
Ideal for use with sulfides, heavy metals, organics, Tris buffer, or protein solids
Simple to calibrate with automatic standard recognition—USA and NIST
Valox® body resists most chemicals and is rugged, dustproof, waterproof,
and IP67 rated

Part Number: WD-35634-35

List: \$36.75

Sale Price: \$31.24

pH calibration buffer pouches assorted
pH 4.00, 7.00 & 10.00 20/pk

Part Number: WD-35653-04



Promo Code: Tester50

Offer expires: July 31, 2019

For more information on these and other fine products, please visit:
http://www.foxscientific.com/pdf/FOX_Oakton.pdf

