


# Additional Shamrock Labels for Cryogenics

## CRYO-STICK LABELS IN COLORS ON SHEETS OR ROLLS -80°C



FCS-43

8 Standard Colors · 5 Fluorescent Colors

Glove Friendly · Smudge Resistant · Supplied in Lots of 500 or 1,000 Labels



FCS-43-4  
(Red)



FCS-43-2  
(Yellow)




FCS-43-30  
(Fl. Chartreuse)

7/16" Dia. Dots (32 Dots per 2 1/2" x 4" Sheet)



FCS-5125-11  
(Rose)




FCS-5125-6  
(Blue)




FCS-5125-3  
(Green)

1/2" x 1 1/4" Labels (32 Labels per 5" x 4" Sheet)



FCS-5125

### STANDARD AND FLUORESCENT COLORS



To order labels in color, please add the color number shown to the label's catalog number.

For instance FCS-43 in fluorescent red would be FCS-43-27. To indicate Packs of 500 add the suffix "D" to the catalog number: FCS-43-27-D

## LABELS FOR ALL CRYOGENIC TEMPERATURES ON SHEETS OR ROLLS


Glove Friendly · Smudge Resistant · Supplied in Lots of 500 or 1,000 Labels

### Cryogenic Labels in Rolls for Thermal Printers



Rolls Cat. No.	Label Size	Materials
KDS-43R FCS-43R LCS-43R	7/16" Dia.	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-59R FCS-59R LCS-59R	1/2" x 9/10"	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-5125R FCS-5125R LCS-5125R	1/2" x 1 1/4"	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-515R FCS-515R LCS-515R	1/2" x 1 1/2"	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-759R FCS-759R LCS-759R	3/4" x 9/10"	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-7512R FCS-7512R LCS-7512R	3/4" x 1 1/4"	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-7515R FCS-7515R LCS-7515R	3/4" x 1 1/2"	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C

### Cryogenic Labels in Sheets for Hand Marking



Sheets Cat. No.	Label Size	Sheet Size	Materials
KDS-43 FCS-43 LCS-43	7/16" Dia.	2 1/4" x 4" 32 labels/sheet	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-59 FCS-59 LCS-59	1/2" x 9/10"	4 1/2" x 4" 40 labels/sheet	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-5125 FCS-5125 LCS-5125	1/2" x 1 1/4"	5" x 4" 32 labels/sheet	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-515 FCS-515 LCS-515	1/2" x 1 1/2"	4" x 4 1/2" 24 labels/sheet	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-759 FCS-759 LCS-759	3/4" x 9/10"	4 1/2" x 4 1/2" 36 labels/sheet	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-7512 FCS-7512 LCS-7512	3/4" x 1 1/4"	4 1/2" x 5" 24 labels/sheet	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-7515 FCS-7515 LCS-7515	3/4" x 1 1/2"	4 1/2" x 4 1/2" 18 labels/sheet	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C

Use catalog numbers above for lots of 1,000 labels. For 500 labels add the suffix "-D" to the appropriate catalog number (KDS-43R-D, KDS-43-D).

### Wrap-Around Labels


These unique labels for sub-zero use feature a matte white area for thermal transfer printing or marking with CryoPen-26.

**Note:** The extended clear portion of each label covers and protects printed information.

The ACTG-WAPC-1 label features an indexing notch between the printable portion of each label.



ACTG-WAPC-1  
1 5/16" x 1 3/4"




WA-1  
1 5/16" x 2 3/8"



### CryoPen-26

Supplied in packs of 4, ink won't fade or smudge.



### Thermal Transfer Ribbons

Shamrock offers a wide variety of ribbons specifically designed to produce an imprint that remains clear and legible throughout cryogenic procedures. See page L-61 or call Shamrock for complete information.

Shamrock Scientific Specialty Systems · 34 Davis Drive, Bellwood, IL 60104

CALL 1-800-323-0249 · FAX 1-800-248-1907 · EMAIL sales@ShamrockLabels.com · SHOP www.ShamrockLabels.com




# Shamrock Cryogenic Labels

## Laser Printer Labels on 8-1/2" x 11" Sheets


### LIQUI-STICK LABEL LASER SHEETS -196°C

- Withstand Liquid or Vapor Phase Liquid Nitrogen
- Resistant to UV and Gamma Irradiation • 7 Color Choices for Labels Starred
- Smudge Resistant Surface for Laser Printing and CryoPen Marking




LASER-CL12T  
(7 Color Choices)


LASER-CL23T  
(7 Color Choices)




LASER-CL48T  
(7 Color Choices)



LASER-CL60 (White)  
for Tubes & Tube Tops



SLS-LASER-CL32T  
(White)



LASER-TRCL12  
(Transparent)



LASER-TRCL23  
(Transparent)

Catalog Number	Label Size	Description	Labels/Sheet	Sheets/Pack
LASER-CL12T	.5" x .94"	7 Color Choices	212	16
★ LASER-CL23T	.5" x 1.28"	7 Color Choices	176	16
★ LASER-CL48T	.4333" Dia.	7 Color Choices	120	16
★ LASER-CL60	.94" x 1.125" & .44" Dia.	White Only	84	16
SLS-LASER-CL32T	.78" x .94"	White Only	70	16
LASER-TRCL12	.5" x .94"	Transparent	84	16
LASER-TRCL23	.5" x 1.28"	Transparent	70	16


★ Seven Colors Available 1 White 2 Yellow 3 Green 4 Red 5 Orange 6 Blue 7 Pink

To indicate label color when ordering, place the color number after the catalog number above.  
For instance LASER-CL12T in yellow would be ordered as LASER-CL12T-2.

### LASER SMUDGEPROOF KIMDURA® SHEETS -54°C

- Withstand Temperatures as low as -54°C • Supplied on 8½" x 11" Sheets/5,000 Labels per Carton
- Glove friendly, smudge resistant • 13 Color Choices


Catalog Number	Label Size	Labels/Sheet	Sheets/Ctn
LASER-KDS-43	7/16" Dia.	212	24
LASER-KDS-59	1/2" x 9/10"	176	29
LASER-KDS-5125	1/2" x 1 1/4"	120	42
LASER-KDS-515	1/2" x 1 1/2"	100	50
LASER-KDS-759	3/4" x 9/10"	104	48
LASER-KDS-7512	3/4" x 1 1/4"	84	60
LASER-KDS-7515	3/4" x 1 1/2"	70	72



#### STANDARD COLORS



#### FLUORESCENT COLORS



To order labels in color, please add the color number shown to the label's catalog number.  
For instance LASER-KDS-43 in fluorescent red would be LASER-KDS-43-27.

Printing template information is supplied with all Laser Labels.