

Amber Glass Wide Mouth Packers

Finneran Products Certified For Science™



-Amber glass bottles ideal for UV light-sensitive liquid samples

-Choice of closure and liner

-Available in standard, precleaned or precleaned/certified in accordance with recommended E.P.A. protocol

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Precleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

CLASS 3: (Certified) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a "Certificate of Analysis."

An independent laboratory is used to certify the container's cleanliness, assuring impartial results.

Standard

Assembled with White Polypropylene Closure, .015" PTFE Lined

Amber Glass Wide Mouth Packers - Class 1

Cat. No.	Description	Qty
9-219	1/2 oz, 15cc Packer, 28-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-220	1 oz, 30cc Packer, 28-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-221	2 oz, 60cc Packer, 33-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-222	4 oz, 120cc Packer, 38-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-223	8 oz, 250cc Packer, 45-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-224	16 oz, 500cc Packer, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-225	32 oz, 950cc Packer, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-226	80 oz, 2500cc Packer, 70-400mm Thread, White Polypropylene Closure, PTFE Lined	4 Pack

Don't see what you are looking for?
Please contact us. Any combination of products is available.

Amber Glass Wide Mouth Packers

Precleaned

Assembled with White Polypropylene Closure, .015" PTFE Lined

Amber Glass Wide Mouth Packers - Class 2

Cat. No.	Description	Qty
9-219-2	1/2 oz, 15cc Packer, 28-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-220-2	1 oz, 30cc Packer, 28-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-221-2	2oz, 60cc Packer, 33-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-222-2	4 oz, 120cc Packer, 38-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-223-2	8 oz, 250cc Packer, 45-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-224-2	16 oz, 500cc Packer, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-225-2	32 oz, 950cc Packer, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-226-2	80 oz, 2500cc Packer, 70-400mm Thread, White Polypropylene Closure, PTFE Lined	4 Pack

Precleaned & Certified

Assembled with White Polypropylene Closure, .015" PTFE Lined

Amber Glass Wide Mouth Packers - Class 3

Cat. No.	Description	Qty
9-220-3	1 oz, 30cc Packer, 28-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-221-3	2 oz, 60cc Packer, 33-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-222-3	4 oz, 120cc Packer, 38-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-223-3	8 oz, 250cc Packer, 45-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-224-3	16 oz, 500cc Packer, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-225-3	32 oz, 950cc Packer, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Amber Glass Wide Mouth Packers with Assorted Closures

Standard

Assembled with Black Phenolic Closure/PE Cone Lined

Amber Glass Wide Mouth Packers - Class 1

Cat. No.	Description	Qty
D0173-1/2	1/2 oz, 15cc Packer, 28-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0173-1	1 oz, 30cc Packer, 28-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0173-2	2 oz, 60cc Packer, 33-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0173-4	4 oz, 120cc Packer, 38-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack

Assembled with White Polypropylene Closure, F217 Lined

Amber Glass Wide Mouth Packers - Class 1

Cat. No.	Description	Qty
D0180-1/2	1/2 oz, 15cc Packer, 28-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0180-1	1 oz, 30cc Packer, 28-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0180-2	2 oz, 60cc Packer, 33-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0180-4	4 oz, 120cc Packer, 38-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0180-8	8 oz, 250cc Packer, 45-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0180-16	16 oz, 500cc Packer, 53-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
D0180-32	32 oz, 950cc Packer, 53-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
D0180-80	80 oz, 2500cc Packer, 70-400mm Thread, White Polypropylene Closure, F217 Lined	4 Pack

Assembled White Polypropylene Closure, Unlined

Amber Glass Wide Mouth Packers - Class 1

Cat. No.	Description	Qty
D0178-1/2	1/2 oz, 15cc Packer, 28-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0178-1	1 oz, 30cc Packer, 28-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0178-2	2 oz, 60cc Packer, 33-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0178-4	4 oz, 120cc Packer, 38-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0178-8	8 oz, 250cc Packer, 45-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0178-16	16 oz, 500cc Packer, 53-400mm Thread, White Polypropylene Closure, Unlined	12 Pack
D0178-32	32 oz, 950cc Packer, 53-400mm Thread, White Polypropylene Closure, Unlined	12 Pack
D0178-80	80 oz, 2500cc Packer, 70-400mm Thread, White Polypropylene Closure, Unlined	4 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Amber Glass Wide Mouth Packers with Assorted Closures

Standard

Without Closures

Amber Glass Wide Mouth Packers - Class 1		
Cat. No.	Description	Qty
D0179-1/2	1/2 oz, 15cc Packer, 28-400mm Thread	24 Pack
D0179B-1/2	1/2 oz, 15cc Packer, 28-400mm Thread - Bulk Pack	624 Pieces
D0179-1	1 oz, 30cc Packer, 28-400mm Thread	24 Pack
D0179B-1	1 oz, 30cc Packer, 28-400mm Thread - Bulk Pack	432 Pieces
D0179-2	2 oz, 60cc Packer, 33-400mm Thread	24 Pack
D0179-4	4 oz, 120cc Packer, 38-400mm Thread	24 Pack
D0179-8	8 oz, 250cc Packer, 45-400mm Thread	24 Pack
D0179-16	16 oz, 500cc Packer, 53-400mm Thread	12 Pack
D0179-32	32 oz, 950cc Packer, 53-400mm Thread	12 Pack
D0179-80	80 oz, 2500cc Packer, 70-400mm Thread	4 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*