



Scientific, Inc.

Save up to 35% with CHEMetrics® COD Vials

Phone & Fax Numbers:

(800) 369-5524 Toll Free

(817) 783-7500 Local

(817) 783-3571 Fax sales@foxscientific.com



No Difference—CHEMetrics® COD Vials can be used directly with any Hach® factory-programmed instrument calibrations, at savings of up to 35%! Request the Performance Data Report showing CHEMetrics COD Vial performance in Hach instrumentation.

- ✓ **USEPA-accepted for wastewater analysis**
- ✓ **Lower cost per test—save up to 35%**
- ✓ **Mercury-free method available**
- ✓ **Compact and eco-friendly packaging**
- ✓ **No need to modify existing equipment or procedures**
- ✓ **Product support with a personal touch**

Cross Reference to Hach® COD Products

COD Reagent Vials*

CHEMetrics Cat. No.	No. of Tests	Range (ppm)	USEPA-accepted	Hach ¹ Equivalent Cat. No.
K-7350S	25	0-150	Yes	21258-25
K-7360S	25	0-1500	Yes	21259-25
K-7370S	25	0-15,000	No	24159-25
K-7355	150	0-150	Yes	21258-15
K-7365	150	0-1500	Yes	21259-15
K-7375	98	0-15,000	No	24159-15

Mercury-free COD Reagent Vials

CHEMetrics Cat. No.	No. of Tests	Range (ppm)	USEPA-accepted	Hach ¹ Equivalent Cat. No.
K-7351S	25	0-150	No	25650-25
K-7361S	25	0-1500	No	25651-25
K-7371S	25	0-15,000	No	28343-25
K-7356	150	0-150	No	N/A
K-7366	150	0-1500	No	25651-15
K-7376	98	0-15,000	No	N/A

*USEPA-accepted COD Vials can be used for NPDES reporting.

¹NOTE: No endorsement by Hach Company is implied or intended.