



# RNT<sup>®</sup>

## Clear Flexible PVC Tubing

**Have everything except your tubing?**



**Let us help!**



- Clear material, with a dense, smooth bore
- Resistant to oxidation, chemicals and bacteria
- Chemical or thermal bonding
- Can be gas or radiation sterilized
- Extreme flexibility permits easy set-ups
- Available in transparent or opaque colors
- Flushes clean with most conventional cleaners and sanitizers
- Can be used with a variety of fittings: radiator clamp or barbed type

<b>Part#:</b>	<b>ID:</b>	<b>OD:</b>	<b>Wall:</b>	<b>Size:</b>	<b>PSI:</b>
T250-000	3/8"	9/16"	3/32"	13	50MAX

**50ft/pks**

**Other sizes available!**

**No high volume minimums!**

**Ask us about volume pricing!**

**We cross reference part numbers!**

• RNT conforms to critical industry standards, having been approved by FDA, USDA, USP Class VI, NSF, 3A and UL. RNT's purity, performance and reliability make it the material of choice when handling liquids, gases and semi-liquids.

8221 East FM 917 • Alvarado, TX 76009

Phone: (800)369-5524 • Fax: (817)783-3571

[sales@foxscientific.com](mailto:sales@foxscientific.com) • [www.foxscientific.com](http://www.foxscientific.com)