

Z446-K Refrigerated Universal Centrifuge



- Max. swing-out capacity:
28 x 50ml, 68 x 15ml, and 20 x microplates
- Temp Range -20°C to +40°C
- Bundle options for tubes or plates!
- Truly universal, 14 rotor options

As low as \$13,500



	Reg.	Bundle
Z446-K-BND	\$15,001	\$14,000
Z446-K-MPBND	\$15,045	\$13,500

BUNDLE ORDERING INFORMATION*

Z446-K-BND Z446-K w/ Swing out rotor for 15/50ml
Z446-K-MPBND Z446-K w/ Swing out rotor for microplates

Z216-MK Refrigerated High Speed Microcentrifuge



As low as \$5,900

- Temp Range -20°C to +40°C
- High speed, up to 21,631xg
- 2 innovative rotor bundle options!



BUNDLE ORDERING INFORMATION*

Z216-K-BND44	Z216-MK Refrigerated w/ 44 place rotor
Z216-K-CMB	Z216-MK Refrigerated w/ COMBI-Rotor™

Reg.	Bundle
\$7,012	\$5,900
\$7,714	\$6,400

Ambient versions also available

Also available:

MC-24™ Touch

High Speed Microcentrifuge

- Combi-rotor for tubes and strips
24 x 1.5/2.0ml + 2 x PCR Strips
- Up to 16,800 xg!

\$2,375

*Additional bundle models available upon request.
for 230V voltage, add "-E"

Offers valid until Dec. 31st 2022

BucketBuddy™

Magnetic Bucket and Tube Rack

Magnetic Modular

Accepts:
Hermle inserts/buckets
up to 750ml
(150 mm x 145 mm area)

13 mm dia.
17 mm dia.
30 mm dia.

Part# Z326-BKR

FREE With ANY Hermle Universal Bundles
(Z306, Z326, Z366, Z446, Z32-HK and Z36-HK Series)

2022 Hermle Centrifuge Bundles



Fox Scientific, Inc.



Unbeatable value!
\$4,500

Compact, only 14in. wide

Z287-A Micro Bundle

- Bundle includes 2 rotors
18 x 1.5 / 2.0 ml and 2 x microplate
- Up to 14,000 rpm and 16,058 xg
- Indicator Light, easily view cycle status
- 11 rotor options!

Bundle Includes



2 x microplate
3,800 rpm / 1,453 xg



18 x 1.5 / 2.0 ml
14,000 rpm / 15,994 xg

BUNDLE ORDERING INFORMATION*

Z287-BND-MB Z287-A with 2-place MP rotor and 18 x 1.5/2.0ml rotor

Reg.	Bundle
\$5,254	\$4,500



As low as \$14,000

Z36-HK SUPER Speed Centrifuge

- Speeds up to 65,395 xg!
- Max. swing-out capacity: 12 x 50ml or 28 x 15ml
- 3 bundle options to meet any need!
- Refrigerated, -20°C to 40°C
- Includes FREE **BenchPACK!™** lab starter pack! (**\$1,500 Value!**)

Bundle Choices:



6 x 250ml rotor
10,000 rpm / 15,650 xg



6 x 50ml rotor
21,000 rpm / 41,410 xg



12 x 1.5/2.0 ml rotor
30,000 rpm / 65,395 xg

BUNDLE ORDERING INFORMATION*

	Reg.	Bundle
Z36-BND250	\$19,718	\$15,800
Z36-BND50	\$17,863	\$14,600
Z36-BND20	\$16,967	\$14,000

Offers valid until Dec. 31st 2022

Additional Bundles Available

Z496-K
Floor Standing*
w/ swing out rotor

\$26,000



*230V Only

SIEVA III
Filtration Centrifuge
w/ perforated or non perforated basket

As low as \$22,000



Z206-A
Compact Centrifuge
w/ 6 x 50ml rotor

\$2,600



Z306
Compact Universal
w/ swing out or microplate options

\$5,400



*Additional bundles available upon request.
To order 230V voltage, add "-E" to part number.