

# Metal Free Centrifuge Tubes

- Validated to contain less than 1ppb of more than 30 trace metal elements, ensuring samples are free of contaminants
- Large white marking area provides excellent contrast for identification
- Graduations are printed in bold, easy-to-read black ink and are accurate to +/-2%
- Easy to write on flat surface polyethylene (HDPE) caps feature integrated leak-proof seal
- Tubes include graduations molded into the conical bottom, ideal for pellet viewing and small volume quantification
- Convenient zip-closure bag is a practical alternative to folding and taping opened bags
- Easily refill reusable racks with tubes packed in bags



**Metal Free Tubes are great for analytical testing, such as...**



**Water**



**Environment**



**Food/Beverage**

Part No.	Description	Qty/Case
<b>230480</b>	15mL Centrifuge Tube, Metal Free, Bag, Sterile	500
<b>230482</b>	50mL Centrifuge Tube, Metal Free, Bag, Sterile	500

Distributed by:



800.369-5524  
 sales@foxscientific.com  
 www.foxscientific.com

