



Chemicals • Laboratory Supplies • Equipment

Scientific, Inc. est. 1988

Ultraseal™ LITE

Semi-Automatic Plate Sealer
and accessories!



Part

Number: Description:

250-10071	Ultraseal LITE Adaptable Semi-Automatic Microplate Sealer 110v/230v
150-10021029	Plate adaptors set of 3
101028121073	Plate adaptor 384 well PCR plates
15010035-102	Spare plate weight to hold film down during sealing - for RECESSED plates only
15010035-101	Spare plate weight to hold film down during sealing - for STANDARD & RIMLESS plates only

Films

Part

Number: Description:

SF01078LS	Transparent Peelable Film 100 sheets 125mm x 78mm
SF01078TS	Transparent Peelable Film - Trial Roll 5 sheets 125mm x 78mm
SF03078LS	Transparent Pierceable Film for Abi 3730 100 sheets 125mm x 78mm
SF03078TS	Transparent Pierceable Film for Abi 3730 - Trial Roll 5 sheets 125mm x 78mm
SF02078LS	Optically Clear Strong Bond Heat Seal for Water Bath PCR 100 sheets 125mm x 78mm
SF02078TS	Optically Clear Strong Bond Heat Seal for Water Bath PCR - Trial Roll 5 sheets 125mm x 78mm
SF04078LS	Peelable Foil (White Laminate) Film for PP Plates 100 sheets 125mm x 78mm
SF04078TS	Peelable Foil (White Laminate) Film for PP Plates - Trial Roll 5 sheets 125mm x 78mm
SF14078LS	Peelable Foil (White Laminate) Film for All Plate Types 100 sheets 125mm x 78mm
SF14078TS	Peelable Foil (White Laminate) Film for All Plate Types - Trial Roll 5 sheets 125mm x 78mm
SF05078LS	Peelable Foil Film for PP Plates, DMSO Resistant 100 sheets 125mm x 78mm
SF05078TS	Peelable Foil Film for PP Plates, DMSO Resistant - Trial Roll 5 sheets 125mm x 78mm
SF06078LS	Easily Pierced Foil Heat Sealing Film 100 sheets 125mm x 78mm
SF06078TS	Easily Pierced Foil Heat Sealing Film - Trial Roll 5 sheets 125mm x 78mm
SF15078LS	Pierceable and Peelable Foil Heat Seal for all Plate Types 100 sheets 125mm x 78mm
SF15078TS	Pierceable and Peelable Foil Heat Seal for all Plate Types - Trial Roll 5 sheets 125mm x 78mm
SF07078LS	Pierceable, Peelable Foil Film to PS Plates 100 sheets 125mm x 78mm

SF07078TS	Pierceable, Peelable Foil Film to PS Plates - Trial Roll 5 sheets 125mm x 78mm
SF08078LS	Strong Bonding Peelable Foil Film for PP Plates 100 sheets 125mm x 78mm
SF08078TS	Strong Bonding Peelable Foil Film for PP Plates - Trial Roll 5 sheets 125mm x 78mm
SF10078LS	Non-woven Gas Permeable Peelable Film for Cell & Seed Culture 100 sheets 125mm x 78mm
SF10078TS	Non-woven Gas Permeable Peelable Film for Cell & Seed Culture - Trial Roll 5 sheets 125mm x 78mm
SF13078LS	60gsm Paper with a Grid Lacquer Coating, Breathable, Peelable Film 100 sheets 125mm x 78mm
SF13078TS	60gsm Paper with a Grid Lacquer Coating, Breathable, Peelable Film - Trial Roll 5 sheets 125mm x 78mm
SF11078LS	Clear Perforated Gas Permeable Peelable Film for Cell & Seed Culture 100 sheets 125mm x 78mm
SF11078TS	Clear Perforated Gas Permeable Peelable Film for Cell & Seed Culture - Trial Roll 5 sheets 125mm x 78mm
229571	Peelable 70um EVA/aluminum laminate sealing foil- color coded red 100 sheets 125mm x 78mm
229572	Easily pierced 20um PP/PS lacquered aluminum foil - color coded green 100 sheets 125mm x 78mm
229573	Pierceable/Peelable 38um lacquered aluminum foil - color coded blue 100 sheets 125mm x 78mm
229574	Strong 85um thick aluminum foil for long term storage - color coded black 100 sheets 125mm x 78mm

8221 East FM 917 • Alvarado, TX 76009
Phone: (800)369-5524 • Fax: (817)783-3571
www.foxscientific.com • sales@foxscientific.com