

Kimberly-Clark\* Exam Gloves provide protection from known and unknown risks while delivering comfort, dexterity and the assurance to do your best science. They exceed industry standards and work hand-in-hand with other head-to-toe solutions from Kimberly-Clark Professional\*. Choose from a wide range of offerings including a true original: our Purple Nitrile\* Exam Gloves, the most recognizable gloves in the market. After use, give your gloves new life through The RightCycle Program, the first large-scale recycling program for non-hazardous glove waste.



	Purple Nitrile*	Sterling* Nitrile	Blue Nitrile	Spring Green Nitrile	Forest Green Nitrile	Lavender* Nitrile
Chemical Protection	<b>HIGHEST</b> Chemotherapy & High Risk Applications		Moderate Chemical Splash		<b>LOWEST</b> Low Risk – Light Splash	
Biohazard Protection	<b>HIGHEST</b> Virus		Blood & Bacteria		<b>LOWEST</b> Fungi	
Tactile Sensitivity	<b>LOWEST</b> Low Sensitivity			<b>HIGHEST</b> High Sensitivity		
Cost	<b>HIGHEST</b>					<b>LOWEST</b>
Food Handling Certified		✓	✓	✓		

### Product Features:<sup>1</sup>

- FDA cleared as examination gloves
- Powder-free
- Not made with natural rubber latex
- Beaded cuff for donning
- Ambidextrous
- Recyclable via The RightCycle Program

<sup>1</sup> Common to All Kimberly-Clark\* Exam Gloves

RightCycle by Kimberly-Clark Professional\* is an innovative recycling program that diverts waste from landfills and helps you achieve, or even exceed, your organization's sustainability goals. The RightCycle Program will turn your used garments and gloves into raw materials to create durable consumer goods. Take the next step in your waste reduction journey and choose RightCycle by Kimberly-Clark Professional\*.



For more information, ask your Distributor Sales Professional or visit [www.kimtech.com/rightcycle](http://www.kimtech.com/rightcycle).

 Products eligible for recycling	K-C Code	Color	Size	Length	Food Handling Certified	Mil Thickness (Fingertip)	AQL	Textured Area	Double Donnable	Total Gloves or Pairs/Case	Gloves or Pairs/Boxes
 Purple Nitrile* Sterile Pairs	55091-55093		S-L	9.5"		5.9	1.0	Fingertips		200 pairs	50 pairs/4
 Purple Nitrile* Sterile Singles	52101-52103		S-L	9.5"		5.9	1.0	Fingertips		400	100/4
 Purple Nitrile-Xtra*	50601-50604; XS = 55090		XS-XL	12"		5.9	1.0	Fingertips		500	50/10
 Purple Nitrile*	55080-55084		XS-XL	9.5"		5.9	1.0	Fingertips		1,000 (XL=900)	100/10 (XL=90/10)
 Smooth Blue Nitrile	50574-50582		XS-XL	9.5"		5.9	1.5			1,000 (XL=900)	100/10 (XL=90/10)
 Textured Blue Nitrile	53100-53104		XS-XL	9.5"		5.9	1.5	Fingertips		1,000 (XL=900)	100/10 (XL=90/10)
 Spring Green Nitrile	43437-43441		XS-XL	9.5"		4.7	0.65	Fingertips		2,000 (XL=1,800)	200/10 (XL=180/10)
 Sterling* Nitrile Sterile Pairs	53130-53132		S-L	9.5"		3.5	1.5	Fingertips		200 pairs	50 pairs/4
 Sterling* Nitrile-Xtra*	53137-53141		XS-XL	12"		3.5	1.5	Fingertips		1,000	100/10
 Sterling* Nitrile	50705-50709		XS-XL	9.5"		3.5	1.5	Fingertips		2,000 (XL=1,700)	200/10 (XL=170/10)
 Lavender* Nitrile	52816-52820		XS-XL	9.5"		2.8	1.5	Fingertips		2,500 (XL=2,300)	250/10 (XL=230/10)
 Forest Green Nitrile	43443-43448		XS-XXL	9.5"		3.5	0.65	Fingertips		2,000 (XXL=1,800)	200/10 (XL=180/10)
PFE-Xtra* Latex <sup>†</sup>	50501-50504		S-XL	12"		10.2	1.5	Fingertips		500	50/10
PFE Latex <sup>†</sup>	57110-57550		XS-XL	9.5"		6.7	1.5	Fingertips		1,000 (XL=900)	100/10 (XL=90/10)
Comfort Latex <sup>†</sup>	43431-43436		XS-XXL	9.5"		5.0	1.5	Fingertips		1,000 (XXL=900)	100/10 (XXL=90/10)