


Gobble up these amazing BLACK FRIDAY DEALS!

**BLACK
FRIDAY
SALE**
THREE WEEKS
LONG

Shake things up!

BEADBUG™
MICROTUBE HOMOGENIZER

- * Powerful mixing action for lysis, grinding and homogenization
- * Eliminates cross contamination
- * Personal sized - 3 tube capacity
- * Cold room safe
- * Wide variety of grinding media available



**BOGO
\$845
VALUE**

D1030	BeadBug 3, 3x2ml tube capacity	\$845.00/2
	Prefilled 2ml tubes	
D1032-30	Triple-Pure™ Zirconium Beads, 3mm, pk 50 tubes	\$158.00
D1032-15	Triple-Pure Zirconium Beads, 1.5mm, pk of 50 tubes	\$158.00
D1032-01	Triple-Pure Zirconium Beads, 0.1mm, pk of 50 tubes	\$158.00
	Additional bead materials/sizes available	

BUY ONE AT \$845 - GET ONE FREE!

Why waste time changing rotors?

myFUGE™ 12

- * COMBI-Rotor™ for tubes and strips
 - 12x1.5/2.0ml tubes
 - 32x0.2ml PCR tubes
 - 4xPCR strips (8x0.2ml)
- * Twice the capacity of most mini centrifuges
- * Near silent operation
- * Starts and stops with closing and opening the lid

C1012	myFuge 12 Mini Centrifuge with COMBI-Rotor	\$390.00/2
C1008-A5-6	Adapter pack, 0.5ml, pack of 6	\$21.50

COMBI-Rotor

BUY ONE AT \$390.00 - GET ONE FREE!

BenchBlotter™
Benchmark Platform Rocker

- * Perfect for blot washing and gel staining
- * Large 12x12 inch platform
- * Gentle mixing - 12rpm, ±8° tilt
- * Cold room and incubator safe
- * Compact, low profile design

BR1000
BR1000-DIMPLED
BR1000-STACK

BenchBlotter Platform Rocker with non-slip mat
Optional dimpled mat for tubes
Stacking platform

\$485.00/2
\$ 38.75
\$111.25

BUY ONE AT \$485.00 - GET ONE FREE!

Call Us Today!
BLACK FRIDAY
DEALS END
12/3/2021




Black Friday BLOWOUT Deals

MY GEL™ InstaView

Run Gels - View Progress - Document Results: All in One Place, All in Real Time


- * Designed for electrophoresis of DNA
- * Blue light illumination - will not damage DNA
- * Capture high resolution images without removing the gel using a smart phone
- * Complete system includes gel box with integrated power supply, casting set, blue light transilluminator and imaging enclosure



E1201 myGel InstaView Electrophoresis System \$875.00
 Choose 2 below, mix or match, with E1201 purchase - a \$216 value
E4500-LD SmartGlow™ Loading Dye with green stain ~~\$108.00~~ FREE!
E4500-PS SmartGlow PreStain with green stain ~~\$108.00~~ FREE!



2 FREE VIALS OF SMARTGLOW WITH PURCHASE - \$216 VALUE!



QuadCount™ Automated Cell Counter

- * User friendly, touch-screen interface
- * Auto-focus with motorized XYZ stage
- * LED light source, precision optics and CMOS detector
- * Quick, normal and precise counting modes
- * Process 4 samples at once with QuadSlide™ disposable slides



Just load the slide and press start


E7500 QuadCount Automated Cell Counter \$3695.00
E7500-S1 QuadSlides, 4 chambered disposable slide, pk of 50 ~~\$89.50~~ 2 FREE PACKS!

STOP WASTING TIME WITH MANUAL COUNTING!

OS-20L Overhead Stirrers

- * Precise digital control, LCD display
- * Feed-through chuck for height adjustment
- * High torque for stirring low to mid viscosity solutions
- * 20 liter stirring capacity
- * Supplied complete with 4 arm propeller and stand assembly

HURRY! SAVE \$160!



IPS2050-20 OS-20L Overhead Stirrer, 20L, with stand assembly	\$859.00	\$699.00
Accessory paddles	Stainless Steel	Teflon
One line propeller	S2050-P-S1 \$89.00	S2050-P-T1 \$139.00
Centrifugal impeller	S2050-P-S2 \$89.00	S2050-P-T2 \$139.00
3 hole paddle	S2050-P-S3 \$89.00	S2050-P-T3 \$139.00
4 arm propeller (included)	S2050-P-S4 \$89.00	S2050-P-T4 \$139.00

DISCOUNTED PRICE! SAVE \$160!



More Black Friday BLOWOUT Deals



myTEMP^{Mini} CO₂

Personal Sized CO₂ Incubator

- * Digital temperature and CO₂ control
- * 13x14.5" footprint (20L)
- * Holds flasks and bottles up to 2L
- * Dual beam IR sensor

H2300-HC2	myTemp Mini CO ₂ Incubator	\$3745.00	\$2999.00
H2300-REG	Optional CO ₂ regulator	\$315.00	
H2300-SH	Extra stainless steel shelf	\$99.00	




DISCOUNTED PRICE! SAVE \$746!

myBLOCKTM HL

With Heated Lid

- * Heated lid prevents evaporation and condensation
- * Compact design - fits in the palm of your hand
- * Digital control of block and lid temperature
- * Fast heating - reaches 37°C in 2 minutes

CHOOSE A FREE BLOCK!



BSH200-HL	myBlock HL Dry Bath with Heated Lid*	\$425.00	\$399.00
Choose ONE of the blocks below, FREE			
BSH100-02	Block, 30x0.2ml or 8 PCR strips	\$89.50	FREE
BSH100-05	Block, 24x0.5ml	\$89.50	FREE
BSH100-15HL	Block, 15x1.5ml	\$89.50	FREE

*Blocks sold separately

Perfect for enzymatic reactions, LAMP assays and other isothermal reactions

DISCOUNTED PRICE AND FREE BLOCK - \$115.50 VALUE!

These offers carry weight!

Accuris Dx Analytical Balances

- * 120 or 220g capacity, readable to 0.1mg
- * Waterproof and wear resistant front panel
- * Multiple weighing modes and print capability
- * Internal calibration

Accuris Precision Balances

- * 320g capacity, readable to 0.01g
- * Multiple weighing modes
- * Easy in-lab calibration



W3101A-120	Accuris Dx Analytical Balance, 120g	\$1850.00	
W3101A-220	Accuris Dx Analytical Balance, 220g	\$2150.00	
W3200-320	Accuris Precision Balance, 320g	\$1045.00	FREE

**FREE Precision Balance with Analytical Balance Purchase!
A \$1045 VALUE!**

Check out these amazing prices!



Z216 MK HERMLE

**\$5500 WITH
44 PLACE
ROTOR!**

- * High speed - max 15,000rpm/21,631xg
- * 11 rotor options for tubes 0.2ml to 5ml
- * Powerful refrigeration - RT to 4° in <10 minutes
- * Compact design: 11x21.6x11.4 inches
- * Complete with 44x1.5/2.0ml rotor



Z216-K-BND44

Z216-MK Refrigerated Microcentrifuge
complete with 44x1.5/2.0ml rotor

\$5500.00

Additional rotors and packages available. Call for details.



44 Place Rotor

\$6684 IF PURCHASED SEPARATELY - SAVE \$1184!

Accuris MR9600 Absorbance Microplate Reader


- * Stand alone operation - no computer needed
 - * User-friendly, touchscreen programming
 - * Integrated plate shaker
 - * Supplied with 405, 450, 492 and 630nm filters
- Additional filters sold separately



**BUNDLE
READER AND
WASHER TO
SAVE BIG!**

Accuris MW9600 Microplate Washer

- * Precision needle position control
- * Auto rinse to eliminate cross-contamination
- * Liquid level sensors in bottles



MR9600-BND

Microplate Reader and Washer Bundle \$8795.00

Choose 1 washer manifold, below - FREE - up to \$325 Value!

MW90000-12CM 12 channel washer manifold

~~\$325.00~~ **FREE**

MW9600-8CM 8 channel washer manifold

~~\$295.00~~ **FREE**



MICROPLATE READER-WASHER BUNDLE - SAVE UP TO \$1325!

hybex™ Erlenmeyer Flasks

- * hybex 3.3 borosilicate glass
- * Autoclavable, reusable and breakage resistant
- * Permanent white enamel graduations and marking area
- * Threaded top with PP GL45 screw cap and drip ring
- * Ideal for shaking, media preparation and storage

Volume	Pack of		
250ml	12	B3220-250	\$229.00
500ml	6	B3220-500	\$139.00
1000ml	6	B3220-1000	\$189.00
2000ml	1	B3220-2000	\$69.00



BUY ONE, GET ONE OF EQUAL OR LESSER VALUE FREE!

**CHECK OUT
THIS
BOGO!**

Free offers require purchase of designated items at listed price. Offer expires 12/03/21. Not valid on previous purchases. Special price applies only to item specified. Other capacities or items cannot be substituted. Use promo code TURKEY21.

Fox Scientific Inc. Phone 800-369-5524 Email: sales@foxscientific.com