



Media Formulation & Transfer for Vaccine Development

Tubing Clamps & Connectors

Catalog No.	Description
F18225-0000	Nylon Screw-Clamp Compressors; Tubing up to 6.4mm (1/4") O.D.
F18227-0000	Acetal Mini Clamps - Tubing under 4.7mm (3/16") O.D.
F18229-0000	Polyester Maxi Clamps - Tubing up to 19mm (3/4") O.D.
F18218-0000	Jaw Style Clamps - Tubing up to 1.3cm (1/2") O.D.
F18212-0000	Screw Style Clamps - Tubing up to 14mm (1/2") O.D.
F19604-0000	"T" Shaped Connectors for 6mm (1/4") Tubing; PP
F19612-0000	Wye (Y) Connectors for 3.2mm (1/8") Tubing; PP
F19616-0000	Wye (Y) Connectors for 9.5mm (3/8") Tubing; PP
H19570-0000	22-Piece Connector Assortment for 3.2 to 12.7mm (1/8 to 1/2") Tubing; PP
F19626-0000	28-Piece Plastic Fitting Kit for 3.2 to 12.7mm (1/8 to 1/2") Tubing



More sizes & styles available

Stopcocks

Catalog No.	Description
F30890-0000	Two-Way for 6.4mm to 9.5mm (1/4" to 3/8") I.D. Tubing; 2mm Bore, PTFE
F30895-0000	Three-Way for 6.4mm to 9.5mm (1/4" to 3/8") I.D. Tubing; 4mm Bore, PTFE
F30888-0010	Three-Way for 10mm (3/8") I.D. Tubing; 8mm Bore, PVDF
F30917-0000	Two-Piece for 12.7mm (1/2") I.D. Tubing; PE

Powder Funnels

Catalog No.	Description
F14660-1065	Sterileware® Wide Stem Sterile PP; 49ml (1.7oz); Pack of 36
F14660-1150	Sterileware® Wide Stem Sterile PP; 784ml (26.5oz); Pack of 24
H14660-0065	Wide Stem PP 49ml (1.7oz); Pack of 12
H14660-0150	Wide Stem PP 784ml (26.5oz); Pack of 4
H14684-0000	Extra Wide Stem 3.8cm (1 1/2"), HDPE 510ml (17.2oz)
H14716-0000	Industrial Drum, Heavy Duty PE 18L (4.75 gal)

ProCulture® Orbital Shaker Platform

Catalog No.	Description
F37041-0000	Turns a magnetic stirrer into an orbital shaker

ProCulture® 37°C Warming Plate

Catalog No.	Description
F37015-0000	120V
F37015-1000	220V, CE Marked

Spinbar® PTFE Octagon Magnetic Stirring Bars

Catalog No.	Description
F37110-0078	7/8 x 5/16" (22.2 x 8mm)
F37110-0001	1 x 5/16" (25.4 x 8mm)
F37110-0003	3 x 1/2" (76.2 x 12.7mm)

Spinbar® PTFE Giant Polygon Magnetic Stirring Bars

Catalog No.	Description
F37118-0004	4 1/4 x 1 1/16" (108 x 27mm) without Pivot Ring
F37118-0006	6 1/4 x 1 1/16" (159 x 27mm) without Pivot Ring

Wide-Mouth Mason Jars - Autoclavable

Catalog No.	Description
F10916-0000	4L (1gal); Pack of 2
F10917-0000	8L (2gal)

Precisionware Narrow-Mouth Bottles - Autoclavable

Catalog No.	Description
F10631-0007	500ml (16oz); 28mm Closure; Pack of 12
F10631-0008	1000ml (32oz); 38mm Closure; Pack of 6

Precisionware Wide-Mouth Bottles - Autoclavable

Catalog No.	Description
F10632-0005	125ml (4oz); 38mm Closure; Pack of 12
F10632-0006	250ml (8oz); 45mm Closure; Pack of 12

Carboys - Autoclavable

Catalog No.	Description
F10794-0025	Without Spigot; 10L (2.6gal)
F10794-0050	Without Spigot; 20L (5.3gal)
F11846-0025	With Spigot; 10L (2.6gal)
F11846-0050	With Spigot; 20L (5.3gal)

Drum and Carboy Funnel - Autoclavable

Catalog No.	Description
H14712-0200	2.1L (0.56gal)
H14712-0250	4.3L (1.14gal)
H14712-0350	14.1L (3.73gal)

Heavy Duty Polyethylene Rectangular Tanks

Catalog No.	Description
H34080-0000	Without Faucet; 11 x 11 x 10 in.
H34090-0000	Without Faucet; 21 x 16 x 14 in.
H34092-0000	With Faucet; 11 x 11 x 10 in.
H34096-0000	With Faucet; 21 x 16 x 14 in.

ProCulture® Fernbach Glass Shaker Flasks

Catalog No.	Description
BP-3200-098	2800ml, Beaded Neck, No Baffles
BP-3225-100	2800ml, GL45, Bottom Edge Baffles
BP-3250-100	2800ml, Delong Neck, Bottom Center Baffles

Order all glassware on sp-wilmadlabglass.com

Contacts:

SP Bel-Art

Bel-Art Products, a part of SP
661 Route 23 South
Wayne, New Jersey 07470 USA
U.S.: +1-800-423-5278
Intl: +1-973-694-0500
cservice.belart@spindustries.com
www.sp-belart.com

SP Wilmad-LabGlass

SP Wilmad-LabGlass, an SP brand
1172 NW Boulevard,
Vineland, NJ 08360 USA
U.S.: +1-800-220-5171
Intl: +1-856-697-3000
cs@wilmad-labglass.com
www.sp-wilmadlabglass.com



More sizes & styles available

ProCulture® Delong Glass Shaker Flasks

Catalog No.	Description
BP-3300-102	125ml, Standard Baffles
BP-3160-306	250ml, No Baffles
BP-3160-314	2000ml, No Baffles
BP-3350-116	4000ml, Deep Baffles

Order all glassware on sp-wilmadlabglass.com

ProCulture® Beaded Rim Glass Shaker Flasks

Catalog No.	Description
BP-3360-102	125ml, Side Baffles
BP-3360-108	1000ml, Side Baffles

Order all glassware on sp-wilmadlabglass.com

ProCulture® Glass Spinner Flasks

Catalog No.	Description
BP-3400-100	125ml, 75mm Neck, 32mm Side Arm
BP-3400-106	1L, 100mm Neck, 45mm Side Arm
BP-3400-114	15L, 100mm Neck, 45mm Side Arm
BP-3400-116	36L, 100mm Neck, 45mm Side Arm

Order all glassware on sp-wilmadlabglass.com



© 2021 Bel-Art Products, a part of SP. All rights reserved.
Bel-Art Products is a part of SP (scientificproducts.com), a manufacturer of biopharma processing equipment, life sciences equipment, scientific glassware and laboratory supplies serving pharmaceutical, biotechnology, educational, food and beverage and industrial markets.
Care has been taken to ensure that the information in this catalog is accurate, however it is not guaranteed to be so and may have changed since last publication. Product dimensions may be rounded for presentation. If specifications are critical to your application please contact techsupport12@spindustries.com to confirm any data.
All common law and registered trademarks referenced in this document are owned by Bel-Art Products/SP or SP Industries, Inc. ProCulture® is a registered trademark of SP Industries, Inc. and used with permission by Bel-Art Products/SP
Wilmad-LabGlass, FTS, Hotpack are registered trademarks of SP Industries, Inc.

