

Reliable Tools

for Enology Labs



Faster, Easier Titrations with Class A Precision

Titrette® bottletop burette



- Meets Class A standards for glass burettes!
- Popular for sulfur dioxide (SO₂) analyses, acid-base titrations, solvent analyses and more
- Suitable for NaOH, KOH, KI, H₂SO₄, KMnO₄, Iodine and more
- Eliminates meniscus reading errors
- Models available with direct data transmission to PC



Description	Cat. No.	2021 List Price
10 mL (standard)	4760141	\$1,092.00
25 mL (standard)	4760151	1,092.00
50 mL (standard)	4760161	1,092.00

Comfortable, Convenient Meniscus Control

macro™ Pipette Controller

- Replaces old-fashioned pipette bulbs with simple, one-finger control
- Rapid or sensitive filling; fast or drop-wise dispensing
- Fits most pipettes 0.1-200 mL


Description	Cat. No.	2021 List Price
Gray	26200	\$88.80
Blue	26202	88.80
Magenta	26203	88.80
Green	26201	88.80



The Standard in Bottletop Dispensers

Dispensette® S dispensers

- Rugged: no seals to wear out!
- Suitable for most laboratory reagents
- Simple to clean & calibrate
- Recirculation valve reduces splashing and helps conserve reagents
- Autoclavable to 250 °F (121 °C)
- Adapters for common bottle threads included
- Optional flexible discharge tube simplifies serial dispensing



Dispensette® S Volume, mL	Without recirculation valve		With recirculation valve	
	Cat. No.	2021 List Price	Cat. No.	2021 List Price
Dispensette® S Digital				
0.1-1	4600310	\$534.80	4600311	\$568.00
0.2-2	4600320	534.80	4600321	568.00
0.5-5	4600330	534.80	4600331	568.00
1-10	4600340	534.80	4600341	568.00
2.5-25	4600350	722.20	4600351	755.40
5-50	4600360	732.80	4600361	766.00
Dispensette® S, Analog-adjustable				
0.1-1	4600100	468.80	4600101	490.00
0.2-2	4600120	468.80	4600121	490.00
0.5-5	4600130	468.80	4600131	490.00
1-10	4600140	468.80	4600141	490.00
2.5-25	4600150	656.20	4600151	677.60
5-50	4600160	672.20	4600161	694.40
10-100	4600170	1,047.20	4600171	1,069.60

Fixed volume models available. For strong acids and nonpolar organics, ask for the BRAND® Dispensette® S Organic dispensers.

food_bev_Enology_bt_2021



Lab Rats Trust BrandTech!

BRANDTECH®
SCIENTIFIC, INC.

Toll Free: 888-522-2726

www.brandtech.com

Reliable Tools for Enology Labs

Precise, Slim, Light & Dependable Stepper

HandyStep® S Repeating Pipette



The HandyStep® S provides unprecedented versatility in a mechanical repeating pipette. It offers 59 different volume settings with accuracy, precision and dependability.

- Comfortable, slim design is lightweight, well-balanced and offers easy one-handed volume selection
- Flexible; accepts third party tips – no more being locked into a sole source
- Tip ejection system allows tip removal without user handling

HandyStep® touch and HandyStep® touch S models also available

Description	Cat. No.	2021 List Price	
HandyStep® S Repeating Pipette	705110		\$402.80
HandyStep® touch Repeating Pipette	705200		1,055.00
PD-Tip™ // size	Non-sterile tips	BIO-CERT® Pure*	
0.1 mL (100 tips)	705700	\$125.20	705730 \$171.60
0.5 mL (100 tips)	705702	125.20	705732 171.60
1.0 mL (100 tips)	705704	125.20	705734 171.60
1.25 mL (100 tips)	705706	125.20	705736 171.60
2.5 mL (100 tips)	705708	125.20	705738 171.60
5 mL (100 tips)	705710	125.20	705740 171.60
10.0 mL (100 tips)	705712	125.20	705742 171.60
12.5 mL (100 tips)	705714	125.20	705744 171.60
25 mL (50 tips)	705716	84.80	— —
25 mL (25 tips)	—	42.40	705746 52.60
50 mL (25 tips)	705718	42.00	705748 52.60
Multi-pack (20 tips each 0.5, 1, 1.25, 2.5, 5, 10 & 12.5 mL)	705720	169.60	— —

*BIO-CERT® is certified sterile, and DNA-, RNase-, ATP- and endotoxin-free.

Comfortable, Convenient Meniscus Control

accu-jet® pro pipette controller

- Provides smooth, intuitive meniscus control
- Simple, sensitive soft blow-out
- Contoured grip for comfort
- Mode switch permits either powered or gravity dispensing

Description	Cat. No.	2021 List Price
accu-jet® pro pipette controller		
Dark Blue	26330	\$426.20
Magenta	26331	426.20
Green	26332	426.20
Royal Blue	26333	426.20



Lightweight and Rugged

Transferpette® S Adjustable Volume Single Channel Pipettes



Capacity, µL	Cat. No.	2021 List Price
0.1-1	705868	\$297.00
0.1-2.5	705869	297.00
0.5-10	705870	297.00
2-20	705872	297.00
5-50	705873	297.00
10-100	705874	297.00
20-200	705878	297.00
100-1000	705880	297.00
500-5000 (0.5-5 mL)	705882	297.00
1000-10000 (1-10 mL)	705884	297.00

- Simple, thumb-tip volume adjustment, even with gloves!
- Shortest stroke length among leading pipettes reduces strain
- Slim design reaches bottom of 15 mL conicals
- Comfortably fits most hands
- Accepts tips from most quality manufacturers
- Fixed volume single channels also available

Transferpette® S -8/-12 Multichannel Pipettes



Capacity, µL	2021 List Price		2021 List Price	
	8 Channel	12 Channel	8 Channel	12 Channel
0.5-10	705900	\$760.00	705920	\$940.00
5-50	705906	760.00	705926	940.00
10-100	705908	760.00	705928	940.00
20-200	705910	760.00	705930	940.00
30-300	705912	760.00	705932	940.00

Transferpette® S -8/-12 pipettes supplied with performance certificate, operating manual, shelf mount, reagent reservoir, pack of replacement O-rings. One-year warranty.



Compact Vacuum Pump for Filtration

VACUUBRAND® ME 1 Vacuum Pump



VACUUBRAND® ME 1 is a compact, high-performance, easy-to-use vacuum pump which is ideal for single and multiple filtrations. An optional vent control valve with dial gauge is available for variable adjustment of pumping speed.

- Small footprint conserves valuable bench space
- Whisper quiet performance & low vibration minimize background noise in the lab
- PTFE diaphragms and valves offer good chemical resistance
- High performance 14 lpm pumping speed & 100 mbar ultimate vacuum

Model	Cat. No.	2021 List Price
ME 1 Vacuum Pump	20721003	\$911.00
Vacuum regulation valve with manometer for ME 1	20696842	194.80

Need Lab Vacuum?

Use the [BrandTech® Vacuum Selection Guide](#) to find the best vacuum pump for a wide range of applications! Go to www.brandtech.com.