

# Reliable Tools



## for Toxicology and Pain Management

### Deep-well Plates



Color	Well volume	Well shape	Well bottom	Qty./Pack	Cat. No.	2021 List Price
<b>Deep-well plates</b>						
96-well, polypropylene	0.3 mL	round	U-bottom	100	701330	\$246.80
96-well, polypropylene	0.5 mL	round	U-bottom	48	701346	377.60
96-well, polypropylene	1.1 mL	round	U-bottom	24	701350	144.20
96-well, polystyrene	1.1 mL	round	U-bottom	32	701352	181.20
96-well, polypropylene, Low Profile	1.2 mL	round	U-bottom	50	701340	287.40
96-well, polypropylene	2.2 mL	square	U-bottom	24	701354	207.80
96-well, polypropylene	2.2 mL	square	V-bottom	50	DW2200	418.00
96-well, polypropylene, Sterile	2.2 mL	square	V-bottom	50	DW2250	534.00
384-well, polypropylene	0.3 mL	square	V-bottom	48	701355	617.00

### Secure Plate Sealing

#### Mat Covers

Description	Autoclavable (121°C for 15 min.)	Qty./Pack	Cat. No.	2021 List Price
Pierceable Sealing Mat for 96-well		50	M2200	\$407.00
Pierceable Sealing Mat for 96-well, Sterile		50	M2250	450.00
Silicone cover for 384 well, deep well plate		50	701357	598.00
Mat for 0.5 mL plates, PP	Yes	50	701358	396.00
Mat for 1.1 mL plates, mod. PE	Yes	24	701360	115.60
Mat for 1.2 mL plates, TPE	Yes	50	701368	312.00
Mat for 2.2 mL plates, EVA	No	24	701362	92.40



#### Sealing Films\*

Application	Material	Qty./Pack	Cat. No.	2021 List Price
Real-Time PCR and Protein Crystallization	Polyolefin, non-tacky	100 sheets	781391ES	\$263.00
Pierceable and Resealable	Vinyl	100 sheets	701396ES	137.80
Sealing film, aluminum	Aluminum	1 roll of 100	781381	231.00
Cold Storage and Automation	Aluminum	100 sheets	701395ES	122.00
DMSO Sample Storage	Polypropylene	100 sheets	701367	335.20
Sealing Paddle for Sealing Films		2 paddles	701381	19.80



\*More available, visit [www.brandtech.com](http://www.brandtech.com)

### Precise, Slim, Light, & Dependable Repeating Pipette

#### HandyStep® S Repeating Pipette



The HandyStep® S provides unprecedented versatility in a mechanical repeating pipette. It offers 59 different volume settings with accuracy, precision and dependability.

- Comfortable, slim design is lightweight, well-balanced and offers easy one-handed volume selection
- Flexible; accepts third party tips – no more being locked into a sole source
- Tip ejection system allows tip removal without user handling

HandyStep® touch and HandyStep® touch S models also available

Description	Cat. No.	2021 List Price	
HandyStep® S Repeating Pipette	705110	\$402.80	
HandyStep® touch Repeating Pipette	705200	1,055.00	
<b>PD-Tip™ // size</b>		<b>Non-sterile tips</b>	<b>BIO-CERT® Pure*</b>
0.1 mL (100 tips)	705700	\$125.20	705730 \$171.60
0.5 mL (100 tips)	705702	125.20	705732 171.60
1.0 mL (100 tips)	705704	125.20	705734 171.60
1.25 mL (100 tips)	705706	125.20	705736 171.60
2.5 mL (100 tips)	705708	125.20	705738 171.60
5 mL (100 tips)	705710	125.20	705740 171.60
10.0 mL (100 tips)	705712	125.20	705742 171.60
12.5 mL (100 tips)	705714	125.20	705744 171.60
25 mL (50 tips)	705716	84.80	— —
25 mL (25 tips)	—	42.40	705746 52.60
50 mL (25 tips)	705718	42.00	705748 52.60
Multi-pack (20 tips each 0.5, 1, 1.25, 2.5, 5, 10 & 12.5 mL)	705720	169.60	— —

\*BIO-CERT® is certified sterile, and DNA-, RNase-, ATP- and endotoxin-free.

# Reliable Tools for Toxicology and Pain Management

## Lightweight and Rugged Pipettes

### Transferpette® S Adjustable Volume Single Channel Pipettes



Capacity, µL	Cat. No.	2021 List Price
0.1-1	705868	\$297.00
0.1-2.5	705869	297.00
0.5-10	705870	297.00
2-20	705872	297.00
5-50	705873	297.00
10-100	705874	297.00
20-200	705878	297.00
100-1000	705880	297.00
500-5000 (0.5-5 mL)	705882	297.00
1000-10000 (1-10 mL)	705884	297.00

- Simple, thumb-tip volume adjustment, even with gloves!
- Shortest stroke length among leading pipettes reduces strain
- Slim design reaches bottom of 15 mL conicals
- Comfortably fits most hands
- Accepts tips from most quality manufacturers
- Fixed volume single channels also available

### Transferpette® S -8/-12 Multichannel Pipettes



Capacity, µL	8 Channel	2021 List Price	12 Channel	2021 List Price
0.5-10	705900	\$760.00	705920	\$940.00
5-50	705906	760.00	705926	940.00
10-100	705908	760.00	705928	940.00
20-200	705910	760.00	705930	940.00
30-300	705912	760.00	705932	940.00

Transferpette® S -8/-12 pipettes supplied with performance certificate, operating manual, shelf mount, reagent reservoir, pack of replacement O-rings. One-year warranty.



## Test Tube Racks



### Polypropylene Test Tube Racks

Sturdy polypropylene racks with alphanumeric positions for easy sample identification. Racks are stackable, non-conductive and will not float in water baths. They are suitable for use in temperatures from -20°C to 90°C. Racks are supplied in packs of five. They are supplied in three pieces for easy and permanent assembly.



#### Test tube racks, PP

Maximum Tube Diameter, mm	Format	Length mm	Width mm	Height mm	White Cat. No.	Blue Cat. No.	Red Cat. No.	2021 List Price
13	6 x 14	265	126	75	4340000	4340001	—	\$82.40
16	5 x 11	265	126	75	4340060	4340061	4340062	82.40
18	5 x 11	265	126	75	4340010	4340011	4340012	82.40
20	4 x 10	265	126	75	4340020	4340021	4340022	82.40
25	4 x 8	265	126	88	4340030	4340031	4340032	82.40
30	3 x 7	265	126	88	4340040	4340041	4340042	82.40

## Certified Volumetric Flasks and Graduated Cylinders

BrandTech® offers BLAUBRAND® USP volumetric instruments that comply with the Class A error limits required by the United States Pharmacopeia.

#### BLAUBRAND® Volumetric Flask, Class A, USP, Certified, Glass Stopper

Volume (mL)	Error limit (+/-mL)	Stopper NS	Qty./Pack	Cat. No.	2021 List Price
5	0.02	10/19	2	36968	\$52.60
10	0.02	10/19	2	36973	52.60
20	0.02	10/19	2	36975	53.00
25	0.03	10/19	2	36977	53.00
50	0.05	12/21	2	36978	56.80
100	0.08	14/23	2	36979	61.40
200	0.1	14/23	2	36980	81.60
250	0.12	14/23	2	36981	85.00
500	0.2	19/26	2	36982	101.00
1000	0.3	24/29	1	36983	81.00
2000	0.5	29/32	1	36984	107.60

#### BLAUBRAND® Graduated Cylinder, Tall Form, Class A, USP, Certified

Volume (mL)	Subdivision (mL)	Error limit (+/-mL)	Height (mm)	Qty./Pack	Cat. No.	2021 List Price
5	0.1	0.05	115	2	<a href="#">32805</a>	\$35.00
10	0.2	0.1	140	2	<a href="#">32808</a>	35.00
25	0.5	0.17	170	2	<a href="#">32820</a>	36.20
50	1	0.25	200	2	<a href="#">32828</a>	41.20
100	1	0.5	260	2	<a href="#">32838</a>	43.00
250	2	1	335	2	<a href="#">32848</a>	82.00
500	5	2.0	365	2	<a href="#">32854</a>	141.40
1000	10	3.0	465	1	<a href="#">32862</a>	99.20
2000	20	6.0	505	1	<a href="#">32864</a>	184.20

For more BLAUBRAND® USP glassware products see [brandtech.com](http://brandtech.com)



Lab Rats Trust BrandTech®!

