

Dispensette® S Bottletop Dispensers

Dispensette S bottletop dispensers build on the fifty year history of BRAND® dispensing expertise. Continual upgrades makes this the safest and most convenient bottletop dispenser ever, all while retaining the features making the Dispensette®

the world's favorite bottletop dispenser. They mount directly on most solvent and reagent bottles for faster, more convenient dispensing. Instruments are autoclavable at 121° C (250° F) for use with sterile reagents.

Easy Calibration Technique

Calibration and adjustments according to ISO 9001 and GLP are done within seconds.

Discharge valve with safety ball closes when discharge tube is not mounted

Discharge tube without recirculation valve

Large viewing window
Easy priming verification

Hinged screw cap
Swings out of the way when dispensing

Calibration mechanism

Volume selection
using interior scalloped track

Discharge tube with recirculation valve

360° rotating valve block with GL 45 thread

Olive-shaped filling valve
for firmer attachment of the filling tube

Recirculation tube
for discharge tubes with recirculation valve

Telescoping filling tube
Easy height adjustment without cutting

Product Features:

Both the Dispensette® S and Dispensette S Organic are constructed using the “floating piston” principle. Each piston is matched individually with precise tolerances to the cylinder of the instrument. A thin film of the dispensed liquid of just a few µm thick acts as a non-wearing seal that reduces friction, so dispensing is easy and convenient.

Dispensette S						Without recirculation valve	With recirculation valve
Volume, mL	Increments, mL	A* < ±		CV* ≤		Cat. No.	Cat. No.
		%	µL	%	µL		
Dispensette S Digital							
0.1-1	0.005	0.6	6	0.2	2	4600310	4600311
0.2-2	0.01	0.5	10	0.1	2	4600320	4600321
0.5-5	0.02	0.5	25	0.1	5	4600330	4600331
1-10	0.05	0.5	50	0.1	10	4600340	4600341
2.5-25	0.1	0.5	125	0.1	25	4600350	4600351
5-50	0.2	0.5	250	0.1	50	4600360	4600361

Dispensette S, Analog-adjustable							
0.1-1	0.02	0.6	6	0.2	2	4600100	4600101
0.2-2	0.05	0.5	10	0.1	2	4600120	4600121
0.5-5	0.1	0.5	25	0.1	5	4600130	4600131
1-10	0.2	0.5	50	0.1	10	4600140	4600141
2.5-25	0.5	0.5	125	0.1	25	4600150	4600151
5-50	1.0	0.5	250	0.1	50	4600160	4600161
10-100	1.0	0.5	500	0.1	100	4600170	4600171

Dispensette S, Fixed-volume							
1		0.6	6	0.2	2	4600210	4600211
2		0.5	10	0.1	2	4600220	4600221
5		0.5	25	0.1	5	4600230	4600231
10		0.5	50	0.1	10	4600240	4600241

Dispensette S Organic						Without recirculation valve	With recirculation valve
Volume, mL	Increments, mL	A* < ±		CV* ≤		Cat. No.	Cat. No.
		%	µL	%	µL		
Dispensette S Organic, Digital							
0.5-5	0.02	0.5	25	0.1	5	4630330	4630331
1-10	0.05	0.5	50	0.1	10	4630340	4630341
2.5-25	0.1	0.5	125	0.1	25	4630350	4630351
5-50	0.2	0.5	250	0.1	50	4630360	4630361

Dispensette S Organic, Analog-adjustable							
0.5-5	0.1	0.5	25	0.1	5	4630130	4630131
1-10	0.2	0.5	50	0.1	10	4630140	4630141
2.5-25	0.5	0.5	125	0.1	25	4630150	4630151
5-50	1.0	0.5	250	0.1	50	4630160	4630161
10-100	1.0	0.5	500	0.1	100	4630170	4630171

Dispensette S Organic, Fixed-volume							
5		0.5	25	0.1	5	4630230	4630231
10		0.5	50	0.1	10	4630240	4630241

Dispensette S Trace Analysis						Without recirculation valve	With recirculation valve
Volume, mL	Valve Spring	A* < ±		CV* ≤		Cat. No.	Cat. No.
		%	µL	%	µL		
Dispensette S Trace Analysis, Analog-adjustable							
1-10	Platinum-iridium	0.5	50	0.1	10	4640040	4640041
1-10	Tantalum	0.5	50	0.1	10	4640240	4640241

A* = Accuracy, CV* = Coefficient of Variation

* The value of accuracy and coefficient of variation are final test values referring to the delivered nominal volume, instrument and distilled water at equilibrium with ambient temperature (20°C) and with smooth operation.



Dispensette S



Dispensette S Organic



Dispensette S Trace Analysis

NEW!



Transferpette® pro

A New Twist on Comfort and Control
Simple. Easy. Efficient

BRANDTECH.com

Transferpette® pro micropipette

The addition of the Transferpette® pro to the trusted Transferpette® family offers scientists even more choice in selecting the ideal pipette for their unique needs. Designed for repeatability and comfort, the Transferpette® pro is built to meet the demands of your daily lab routine.

The Transferpette® pro features a slim, ergonomic handle paired with a rotating finger rest- a unique design that adjusts to fit any hand position, left- or right-handed. Combined with two intuitive calibration systems— Easy Calibration for tool-free, permanent adjustment and User Adjustment for quick, reversible fine-tuning—the Transferpette® pro ensures precise performance, even with changing liquid types, temperatures, or style of tips.

Your benefits:

- + **Simple:** Comfortable pipetting with a slim profile and adjustable finger rest
- + **Easy:** Adapt quickly to changing conditions with new user adjustment
- + **Efficient:** Reliable, durable, and convenient



A New Twist on Comfort and Control



Every application has different requirements. The Transferpette® pro adapts perfectly to your individual needs with additional innovative features:

Adjustable finger rest



Adjustable to 9 positions for an ergonomic hand position

User Adjustment



Easy and temporary adjustment for special conditions such as temperature, tip shape, or liquid type

Counter protection



Protects the pipette by preventing over-rotation, ensuring lasting performance and reliability

Extended volume range



Now available in 300 µl, 2.5 ml, and 1250 µl nominal volume

Counter protection

Prevents damage from over-rotating beyond volume range

Rotatable finger rest

Simple adjustment without tools

User Adjustment

Temporary adjustment to changing conditions

Extended volume range

13 single-channel pipettes from 0.1 µl to 10 ml

6 multi-channel pipettes from 0.5 µl to 1250 µl



All product information can be found at BRANDTECH.com/pro

Transferpette® pro and Transferpette® S in comparison

	Transferpette® pro	Transferpette® S
Volume range single-channel (adjustable)	12 sizes: from 0.1 µl to 10 ml	10 sizes: from 0.1 µl to 10 ml
Volume range multi-channel	6 sizes: from 0.5 µl to 1250 µl	5 sizes: from 0.5 µl to 300 µl
Size of handle	Slim	Standard
Adjustable finger rest	✓	-
One-handed operation	✓	✓
Uses universal tips	✓	✓
Volume lock	✓	✓
Counter protection	✓	-
Easy Calibration Technology	✓	✓
Temporary User Adjustment	✓	-
Completely autoclavable (121 °C)	✓	✓
UV-resistant	✓	✓
Warranty for single-channel pipettes	3 years	3 years
Warranty for multi-channel pipettes	1 year	1 year



(800) 369-5524 • sales@foxscientific.com • www.foxscientific.com

BRANDTECH® Scientific | 888-522-2726 | www.brandtech.com



Product Launch DEAL

NEW!

Transferpette® *pro* Pipettes

Get a **FREE pipette** when you purchase 3 Transferpette *pro* single channel pipettes or 2 Transferpette *pro* multi channel pipettes **up to a \$1,525 USD value!**



Capacity, µL	BRANDTECH Cat. No.	Recommended tip size, µL
Single-Channel Adjustable Pipettes		
0.1 - 1	706868	0.1 - 20
0.1 - 2.5	706869	0.5 - 20
0.5 - 10	706870	0.5 - 20
2 - 20	706871	0.5 - 20
2 - 20	706872	2 - 200
5 - 50	706873	2 - 200
10 - 100	706874	2 - 200
20 - 200	706878	2 - 200
30 - 300	706879	5 - 300
100 - 1,000	706880	50 - 1,000
250 - 2,500	706881	500 - 5,000
500 - 5,000	706882	500 - 5,000
1,000 - 10,000	706884	1,000 - 10,000

Capacity, µL	BRANDTECH Cat. No.	Recommended tip size, µL
8-Channel Adjustable Pipettes		
0.5 - 10	706900	0.5 - 20
5 - 50	706906	2 - 200
10 - 100	706908	2 - 200
20 - 200	706910	2 - 200
30 - 300	706912	5 - 300
50 - 1,250	706916	50 - 1,250
12-Channel Adjustable Pipettes		
0.5 - 10	706920	0.5 - 20
5 - 50	706926	2 - 200
10 - 100	706928	2 - 200
20 - 200	706930	2 - 200
30 - 300	706932	5 - 300
50 - 1,250	706936	50 - 1,250



Redemption Instructions

1. Purchase three Transferpette *pro* single channel pipettes or purchase two Transferpette *pro* multi channel pipettes
2. Complete the online form to select your free item (of equal or lesser value than the items purchased) and upload your proof of purchase at brandtech.com/promoform



All submissions must be completed online and include a valid proof of purchase dated between 10/27/25 - 12/31/26. Allow 4-6 weeks for delivery.

Product appearance, catalog numbers, prices, specifications, and technical information are subject to change without notice. Free item must be of equal or lesser value than any one item purchased. Offer only valid to customers in the US and Canada. Not valid for resellers. Void where prohibited by law or company policy. This offer may not be combined with any other offer. Limits may apply.

888-522-2726 | brandtech.com

2025 Deals

Dispensette® S Bottletop Dispenser

Get a FREE Dispenser when you purchase 3 Dispensette® S Bottletop Dispensers, **up to a \$1,490 USD value!**



Purchase products between **1/1/25 - 12/31/25**. For product information, visit brandtech.com, or call Customer Service at 888-522-2726.

Volume, mL	BRANDTECH Cat. No.	BRANDTECH Cat. No.
	Dispensette® S, Digital	Digital Organic
0.1-1	4600311	—
0.2-2	4600321	—
0.5-5	4600331	4630331
1-10	4600341	4630341
2.5-25	4600351	4630351
5-50	4600361	4630361
	Dispensette® S, Analog	Analog Organic
0.1-1	4600101	—
0.2-2	4600121	—
0.5-5	4600131	4630131
1-10	4600141	4630141
2.5-25	4600151	4630151
5-50	4600161	4630161
10-100	4600171	4630171
	Dispensette® S, Fixed Vol.	Fixed Vol. Organic
1	4600211	—
2	4600221	—
5	4600231	4630231
10	4600241	4630241

Items listed above come with recirculation valve.

For a complete listing of dispensers see brandtech.com



Redemption Instructions

1. Purchase three Dispensette S Bottletop Dispensers
2. Complete the online form to select your free item (of equal or lesser value than the items purchased) and upload your proof of purchase at brandtech.com/promoform



All submissions must be completed online and include a valid proof of purchase dated between 1/1/25 – 12/31/25. Allow 4-6 weeks for delivery.

Product appearance, catalog numbers, prices, specifications, and technical information are subject to change without notice. Offer only valid to customers in the US and Canada. Not valid for resellers. Void where prohibited by law or company policy. This offer may not be combined with any other offer. Limits may apply.

888-522-2726 | brandtech.com