




Laboratory Torches

A properly matched torch can help labs improve day to day reliability, reduce avoidable maintenance, and maintain confidence when working with challenging sample types.



For your ICP needs!



Part #	Size	Description	Instrument Compatability
750-03*PG	1.8mm	Spectroflame Fixed SOP (Side-on-Plasma) Torch	Spectro Modula, CIROS, Genesis, & ARCOS
750-05*PG	1.0mm	Spectroflame Fixed SOP (Side-on-Plasma) Torch	Spectro Modula, CIROS, Genesis, & ARCOS
800-12*PG	Tapered, 2.3mm	Exact match OEM replacement suitable for axial orientations, 90% bend	Agilent / Varian: 700-ES, Liberty, Vista Series
900-51*PG	2.5mm	Recommended torch for aqueous samples, not Hydrofluoric Acid (HF) Compatible	Agilent 7700/7800/7900/8800 /8900 Series ICP-MS